

# Package ‘itraxR’

November 7, 2022

**Type** Package

**Title** Itrax Data Analysis Tools

**Version** 1.8

**Date** 2022-11-07

**Author** Thomas Bishop

**Maintainer** Thomas Bishop <tombishopemail@gmail.com>

**URL** <https://thomasbishop.uk>

**BugReports** <https://github.com/tombishop1/itraxR/issues>

**Description** Parse, trim, join, visualise and analyse data from Itrax sediment core multi-parameter scanners manufactured by Cox Analytical Systems, Sweden. Functions are provided for parsing XRF-peak area files, line-scan optical images, and radiographic images, alongside accompanying metadata. A variety of data wrangling tasks like trimming, joining and reducing XRF-peak area data are simplified. Multivariate methods are implemented with appropriate data transformation.

**License** GPL-3

**Encoding** UTF-8

**RoxygenNote** 7.2.1

**LazyData** true

**Depends** dplyr, ggplot2, R (>= 3.5.0), compositions, grid

**Imports** readr, tiff, janitor, ggcorrplot, rlang, tidyr, broom, tibble, stringr, munsellinterpol

**Suggests** magrittr

**NeedsCompilation** no

**Repository** CRAN

**Date/Publication** 2022-11-07 22:10:02 UTC

## R topics documented:

CD166_19_S1 . . . . .	2
itrax_correlation . . . . .	3
itrax_image . . . . .	4
itrax_import . . . . .	5
itrax_join . . . . .	6
itrax_meta . . . . .	6
itrax_munsell . . . . .	7
itrax_ordination . . . . .	8
itrax_qspecsettings . . . . .	9
itrax_radiograph . . . . .	9
itrax_reduce . . . . .	10
itrax_restspectra . . . . .	11
itrax_section . . . . .	12
itrax_spectra . . . . .	13
<b>Index</b>	<b>15</b>

---

CD166\_19\_S1

*Itrax core scanner data for core CD166\_19\_S1*

---

### Description

A dataset containing XRF, optical and radiographic images and metadata for ocean core CD166\_19\_S1.

### Usage

CD166\_19\_S1

### Format

A list containing the following:

**xrf** a tibble of scan variables

**rgb** a list containing the optical image matrix and associated metadata

**rad** a list containing the radiographic image matrix and associated metadata

@source Wynn, R. B., and B. T. Cronin. 2005. "RRS "Charles Darwin" Cruise CD166, 29 Oct - 22 Nov 2004. Sedimentary processes and deposits in the Agadir Basin and Gulf of Cadiz." 59. Vol. 44.

---

itrax_correlation	<i>Calculate a correlation matrix for Itrax result data</i>
-------------------	---

---

### Description

Calculates a correlation matrix for Itrax data results including normalisation and visualisation

### Usage

```
itrax_correlation(  
  dataframe,  
  elementsonly = TRUE,  
  zeros = "addone",  
  transform = TRUE,  
  plot = FALSE  
)
```

### Arguments

dataframe	pass the name of a dataframe parsed using "itrax_import()" or "itrax_join()"
elementsonly	if TRUE, only chemical elements are included. If FALSE, the data is passed unfiltered, otherwise a character vector of desired variable names can be supplied
zeros	if "addone", adds one to all values. If "limit", replaces zero values with 0.001. Otherwise a function can be supplied to remove zero values.
transform	binary operator that if TRUE will center-log-transform the data, if FALSE will leave the data untransformed. Otherwise, a function can be supplied to transform the data.
plot	set to true if a biplot is required as a side-effect

### Value

a correlation matrix object

### Examples

```
itrax_correlation(CD166_19_S1$xf, plot = TRUE)
```

---

`itrax_image`*Read an Itrax Image File*

---

## Description

Reads an Itrax image file and trims it according to the metadata provided.

## Usage

```
itrax_image(  
  file = "optical.tif",  
  meta = "document.txt",  
  plot = FALSE,  
  trim = TRUE  
)
```

## Arguments

<code>file</code>	defines the name of the datafile to parse
<code>meta</code>	defines the relating metadata
<code>plot</code>	would you like to create a plot as a side-effect?
<code>trim</code>	defines custom trim parameters. The default behaviour uses the limits from the metadata file. Set the false for no trimming, or set the position limits by passing a two element vector.

## Value

a matrix of RGB values, and the relevant data from the metadata file relating to the image.

## Examples

```
itrax_image(file = system.file("extdata",  
                               "CD166_19_S1_optical_lowres.tif",  
                               package = "itraxR",  
                               mustWork = TRUE),  
            meta = system.file("extdata",  
                               "CD166_19_S1_xrf_document.txt",  
                               package = "itraxR",  
                               mustWork = TRUE),  
            plot = TRUE)
```

---

itrax_import	<i>Import Itrax core-scanner result file</i>
--------------	--

---

### Description

Imports and parses data from a results file created by Q-Spec software, part of the Itrax core scanner.

### Usage

```
itrax_import(  
  filename = "Results.txt",  
  depth_top = NA,  
  trim_top = 0,  
  trim_bottom = 0,  
  parameters = "some"  
)
```

### Arguments

filename	defines the name of the datafile to parse
depth_top	defines the coring in depth of the top of the core, in mm
trim_top	defines the length of any trimming required of data at the top of the core, in mm
trim_bottom	defines the length of any trimming required at the bottom of the core, in mm
parameters	one of 'all' (leave all parameters), 'some' (remove some less useful parameters)

### Value

a tibble of the parsed Itrax data

### Examples

```
itrax_import(  
  filename = system.file("extdata",  
    "CD166_19_S1_Results.txt",  
    package = "itraxR",  
    mustWork = TRUE),  
  depth_top = 0)
```

---

itrax_join	<i>Join two or more Itrax result datasets</i>
------------	---

---

**Description**

Join two or more Itrax datasets that have been parsed using "itrax\_import()"

**Usage**

```
itrax_join(list)
```

**Arguments**

**list** a list of dataframes that are parsed Itrax result files — this should have been imported using "itrax\_import()" and must have a depth variable present. This should take the form "list(core1 = core1, core2 = core2)"

**Value**

a tibble of all the input data

**Examples**

```
itrax_join(list(core1 = CD166_19_S1$xrf, core2 = CD166_19_S1$xrf))
```

---

itrax_meta	<i>Parse Itrax scan metadata</i>
------------	----------------------------------

---

**Description**

Parses the "document.txt files" generated from Itrax core scanners

**Usage**

```
itrax_meta(datafile = "document.txt")
```

**Arguments**

**datafile** a "document.txt files" generated from an Itrax core scanner

**Value**

a dataframe of all the parsed input data

**Examples**

```
itrax_meta(system.file("extdata",
                      "CD166_19_S1_xrf_document.txt",
                      package = "itraxR",
                      mustWork = TRUE))
```

---

itrax_munsell	<i>Convert an Itrax Image File into Munsell Colour</i>
---------------	--

---

**Description**

Reads a colour calibrated Itrax image file and processes it to estimate Munsell colour.

**Usage**

```
itrax_munsell(image, proportion = 0.1)
```

**Arguments**

image	defines the name of the image file imported using ‘itrax_image()’. It is essential that the image has been colour calibrated using a colour card or other method.
proportion	defines the width down the centre of the image to use for processing

**Value**

a table of values

**Examples**

```
## Not run:
itrax_image(file = system.file("extdata",
                              "CD166_19_S1_optical_lowres.tif",
                              package = "itraxR",
                              mustWork = TRUE),
            meta = system.file("extdata",
                              "CD166_19_S1_xrf_document.txt",
                              package = "itraxR",
                              mustWork = TRUE),
            plot = FALSE) %>%
magrittr::extract2(1) %>%
itrax_munsell() %>%
dplyr::slice_sample(n = 10)

## End(Not run)
```

---

itrax\_ordination      *Principle Component Analysis on Itrax scan data*

---

### Description

Performs and visualises principle component analysis data from Itrax result data

### Usage

```
itrax_ordination(  
  dataframe,  
  elementsonly = TRUE,  
  zeros = "addone",  
  transform = TRUE,  
  return = "list",  
  plot = FALSE  
)
```

### Arguments

dataframe	pass the name of a dataframe parsed using "itrax_import()" or "itrax_join()"
elementsonly	if TRUE, only chemical elements are included. If FALSE, the data is passed unfiltered, otherwise a character vector of desired variable names can be supplied
zeros	if "addone", adds one to all values. If "limit", replaces zero values with 0.001. Otherwise a function can be supplied to remove zero values.
transform	binary operator that if TRUE will center-log-transform the data, if FALSE will leave the data untransformed. Otherwise, a function can be supplied to transform the data.
return	if "pca" the output of prcomp() is returned, otherwise "list" is a list including the transformed data, sample scores, and the output of prcomp().
plot	set to true if a biplot is required as a side-effect

### Value

either an output of prcomp(), or a list including the input data

### Examples

```
itrax_ordination(CD166_19_S1$xrf, plot = TRUE)
```



---

itrax\_qspecsettings     *Read a Q-Spec settings file and parse the key-value pairs*

---

**Description**

This is used to retrieve settings important elsewhere, for example the mca bin width and offset

**Usage**

```
itrax_qspecsettings(filename = "Results_settings.dfl")
```

**Arguments**

filename            the \*.dfl settings file that relates to the rest of the data

**Value**

a tibble of the parsed data

**Examples**

```
itrax_qspecsettings(filename = system.file("extdata",  
                                           "Results_settings.dfl",  
                                           package = "itraxR",  
                                           mustWork = TRUE)  
)
```

---

itrax\_radiograph     *Read an Itrax Radiograph File*

---

**Description**

Reads an Itrax radiograph file and trims it according to the metadata provided.

**Usage**

```
itrax_radiograph(  
  file = "radiograph.tif",  
  meta = "document.txt",  
  plot = FALSE,  
  trim = TRUE  
)
```

**Arguments**

file	defines the name of the datafile to parse
meta	defines the relating metadata
plot	would you like to create a plot as a side-effect?
trim	defines positions of the trim if required, input is a vector with min and max positions

**Value**

a matrix of RGB values, and the relevant data from the metadata file relating to the image. Also computes the aspect ratio of the image.

**Examples**

```
itrax_radiograph(file = system.file("extdata",
  "CD166_19_S1_radiograph_adj.tif",
  package = "itraxR",
  mustWork = TRUE),
  meta = system.file("extdata",
  "CD166_19_S1_rad_document.txt",
  package = "itraxR",
  mustWork = TRUE),
  plot = TRUE)
```

---

 itrax\_reduce

*Reduce Itrax XRF data*


---

**Description**

Reduces Itrax XRF data into arbitrary chunks using an arbitrary function. This is useful when making direct comparisons between the Itrax XRF data and some other data collected at a lower resolution.

**Usage**

```
itrax_reduce(
  dataframe,
  names = c(1:length(breaks_lower)),
  breaks_lower,
  breaks_upper,
  fun = mean,
  edges = c(">=", "<"),
  by = NULL
)
```

**Arguments**

dataframe	defines the name of the XRF data to reduce, usually a itraxR::itrax_import like tibble
names	optional, a vector of the same length as 'breaks'
breaks_lower	a vector of the lower limit of each chunk
breaks_upper	a vector of the upper limit of each chunk
fun	the function to apply in order to reduce the data. Default is mean(), but sd() is also common
edges	a vector of length 2 with the upper and lower bound behaviour; can be any of '<', '<=', '>', '>='
by	if contiguous samples of even sizes are required, 'by' defines the chunk size and will automatically generate 'breaks'

**Value**

a tibble with the same number of rows as 'breaks' and the same number of columns as 'dataframe'

**Examples**

```
itrax_reduce(dataframe = CD166_19_S1$xf, by = 10)
```

---

itrax\_restspectra      *Make a spectrograph from raw Itrax data spectra files*

---

**Description**

Parses a folder full of raw spectra files from an Itrax core scanner and produces a spectral graph of all the data by position

**Usage**

```
itrax_restspectra(
  foldername = "XRF data",
  parameters = "settings.dfl",
  datapos = 37,
  plot = TRUE,
  trans = "pseudo_log"
)
```

**Arguments**

foldername	defines the folder where the spectra "*.spe" files are located
parameters	optionally, defines the Q-Spec settings file from which to calculate the channel energies
datapos	defines the row at which spectral data begins in the files
plot	TRUE/FALSE, selects whether to create a plot as a side-effect
trans	transformation applied in the plot - see '?ggplot2::scales_colour_gradient()' for options

**Value**

a dataframe of all the spectral data

**Examples**

```
## Not run: itrax_restspectra("~/itraxBook/CD166_19_(2020)/CD166_19_S1/CD166_19_S1/XRF data")
```

---

itrax_section	<i>Cluster analysis and statistical grouping of Itrax data</i>
---------------	--

---

**Description**

Performs a cluster analysis and automatic statistical grouping of parsed Itrax results data to n groups. Also provides information on the most "representative" (central) of each group. These can be used to develop a sub-sampling regime for calibration using another method.

**Usage**

```
itrax_section(
  dataframe,
  divisions = 30,
  elementonly = TRUE,
  zeros = "addone",
  transform = TRUE,
  plot = FALSE
)
```

**Arguments**

dataframe	pass the name of a dataframe parsed using "itrax_import()" or "itrax_join()" or "itrax_reduce()".
divisions	the number of groups to slice into - also the number of representative samples returned.
elementonly	if TRUE, only chemical elements are included. If FALSE, the data is passed unfiltered, otherwise a character vector of desired variable names can be supplied.

zeros	if "addone", adds one to all values. If "limit", replaces zero values with 0.001. Otherwise a function can be supplied to remove zero values.
transform	binary operator that if TRUE will center-log-transform the data, if FALSE will leave the data untransformed. Otherwise, a function can be supplied to transform the data.
plot	set to true if a summary plot is required as a side-effect - the input dataset must have a depth or position variable - depth is used preferentially.

**Value**

the input data with additional columns 'group' and 'calib\_sample', and possibly 'uid' if not supplied.

**Examples**

```
itrax_section(CD166_19_S1$xf, plot = TRUE)
itrax_section(CD166_19_S1$xf %>% itrax_reduce(by = 10), plot = TRUE)
```

---

itrax_spectra	<i>Import an individual spectra file</i>
---------------	--

---

**Description**

Sometimes it is helpful to read an individual spectral file for diagnostics

**Usage**

```
itrax_spectra(filename, parameters = "settings.dfl", plot = TRUE)
```

**Arguments**

filename	defines the name of the *.spe datafile from the core scanner to parse
parameters	optionally defines a relevant Q-Spec settings file in order to compute the energy scale, otherwise channel numbers are reported
plot	logical, if TRUE a side-plot is created

**Value**

a tibble of the parsed data



# Index

## \* datasets

CD166\_19\_S1, [2](#)

CD166\_19\_S1, [2](#)

itrax\_correlation, [3](#)

itrax\_image, [4](#)

itrax\_import, [5](#)

itrax\_join, [6](#)

itrax\_meta, [6](#)

itrax\_munsell, [7](#)

itrax\_ordination, [8](#)

itrax\_qspecsettings, [9](#)

itrax\_radiograph, [9](#)

itrax\_reduce, [10](#)

itrax\_restspectra, [11](#)

itrax\_section, [12](#)

itrax\_spectra, [13](#)